

## CLINICAL STIMULATOR REFERENCE CHART

| MANUFACTURER     | MODEL #                                       | UNIT MADE PRIOR TO 2000    |  | COMPLIANT UNITS                        |
|------------------|---|----------------------------|--|--|
|                  |   | ADAPTOR #                  | WIRE #                                 | WIRE #                                 |
| Amrex            | HV752   | AD4500                     | W4872 (Ban.-Ban.),<br>W4772 (Ban.-Pin) | W4872 (Ban.-Ban.),<br>W4772 (Ban.-Pin) |
|                  | LVG325A                                       | AD4500                     | W4872 (Ban.-Ban.),<br>W4772 (Ban.-Pin) | W4872 (Ban.-Ban.),<br>W4772 (Ban.-Pin) |
|                  | MS Series (all)                               | AD4500                     | W4872 (Ban.-Ban.),<br>W4772 (Ban.-Pin) | W4872 (Ban.-Ban.),<br>W4772 (Ban.-Pin) |
|                  | Z-Stim IF Series                              | AD4500                     | W4872 (Ban.-Ban.),<br>W4772 (Ban.-Pin) | W4872 (Ban.-Ban.),<br>W4772 (Ban.-Pin) |
| Chattanooga      | CPS 200 Combo, CPS 200 Stim                   | No CMP solution for U.S.A. | W18112 (Cdn. Only)                     | W7172                                  |
|                  | CPS 400 Combo, CPS 400 Stim                   | No CMP solution for U.S.A. | W18112 (Cdn. Only)                     | W7172                                  |
|                  | Forte 200, Forte 400 Combo                    | No CMP solution for U.S.A. | W18112 (Cdn. Only)                     | W7172                                  |
|                  | Forte CB 450, Forte ES 450                    | AD8500                     | W6572                                  |  |
|                  | Intelect 120, 170, 775MP, 750MP, 600, 900     | AD8500                     | W6572                                  |  |
|                  | Intelect 500, 500S, 700, 700C, 775            |                            | W2372                                  |  |
|                  | Intelect 550, 551                             | AD4500                     | W4472 or W4372                         |  |
|                  | Intelect Combo, Legend Stim                   | No CMP solution for U.S.A. | W18112 (Cdn. Only)                     | W7172                                  |
|                  | Intelect IFC2 (squared)                       | AD8500                     | W6572                                  |  |
|                  | Intelect Legend XT, Vectra Genisys, Transport |                            |  | W17672                                 |
|                  | Intelect Low Volt 110                         |                            |  | W7972                                  |
|                  | Vectra Pro, Vectra                            |                            |  | W7172                                  |
|                  | VMS   |                            | W2372                                  |  |
|                  | VMS II  | AD8500                     | W6572                                  |  |
| Dynatronics      | Dynatron 50, 125, 525                         | AD6500                     | W6672                                  | W6672                                  |
|                  | Dynatron Plus 550, 650, 850, 950              | AD6500                     | W6672,<br>W6272 (combo)                | W6672,<br>W6272 (combo)                |
| Enraf            | Endomed 182, 482                              |                            |  | W7772 (DIN-pins)<br>W7872 (DIN-ban.)   |
|                  | Endomed 582, 982                              |                            |  | W8372 (DIN-pins)<br>W8272 (DIN-ban.)   |
|                  | Endomed CV-405                                |                            |  | W8596B                                 |
|                  | Sonopulse 464                                 |                            |  | W2872IVBRB                             |
| Enraf (Henley)   | Endomed 433                                   | AD8400                     | W8472 (DIN-pins)<br>W8572 (DIN-ban.)   |  |
| Excel            | Ultra / Multi / Isotron Series                | AD6500                     | W6572,<br>W6372 (combo)                | W6572,<br>W6372 (combo)                |
| Fisher           | Model 700                                     | AD4500                     | W4772 (Ban.-Pin)                       |  |
| H Wave           |   | AD4500                     | W4472 or W4372                         |  |
| Hakomed          |   |                            |  | W8672                                  |
| Healthcare Mnfg. | Masterstim                                    | AD6500                     | W6572                                  | W6572                                  |
| Kaladar (KMI)    | The Stim                                      | AD6500                     | W6572                                  |  |
| LSI              | System II and IV, MS300, MS400                | AD7500                     | W7572                                  | W7272                                  |
| Medco            | Medcollator                                   | AD8500                     | W6572                                  |  |
| Metron           | All models                                    | AD4500                     | W4772 (Ban.-Pin)<br>W4872 (Ban.-Ban.)  | W4772 (Ban.-Pin)<br>W4872 (Ban.-Ban.)  |

## CLINICAL STIMULATOR REFERENCE CHART

| MANUFACTURER              | MODEL #                      | UNIT MADE PRIOR TO 2000       |  | COMPLIANT UNITS                         |
|---------------------------|------------------------------|-------------------------------|--|---|
|                           |                              | ADAPTOR #                     | WIRE #   | WIRE #                                  |
| Mettler                   | Sonicator Plus 992, 994, 930 | AD9500                        | W7472,<br>W7372 (combo)  | W7472,<br>W7372 (combo)                 |
|                           | Sys*Stim 294, 226, 206       | AD9500                        | W7472  | W7472                                   |
| Monad                     | Myomatic                     | AD4500                        | W4472  |   |
| Multiplex                 |                              | AD8400                        | W8072  |   |
| Nemectron                 | Edit                         |                               |  | W8672                                   |
| Neurocare                 | Neurocare 1000               |                               |  | W8972IVBRB                              |
| Omnisound                 | 2500C                        | AD4500                        | W4472BB  |   |
| Orthion Corp.             | Sentry Mark V Model 9000     |                               | W2372  |   |
| P.E.T. Eurotech           | Magnum 100                   | AD6500                        | W6572  |   |
| Phyacton                  |                              |                               |  | W4772 (Ban.-Pin)<br>W4872 (Ban.-Ban.)   |
| Physio Tech.              | Omnistim                     | AD6500                        | W6572  |   |
| Physiodynamic             | Therastim 1000               | AD6500                        | W6412  | W7412                                   |
| Physiomed                 |                              | AD8400                        | W8072  |   |
| Precision<br>Microcurrent |                              | AD4500                        | W4472  | W4472                                   |
| Prostim                   | ET 3000                      |                               |  | W18672                                  |
| Rebox Unit                |                              | AD4500                        | W4472  |   |
| Richmar                   | HV 1000, 2000                | AD8500 with<br>or AD9500 with | W6572<br>W7572   |   |
|                           | HV Combo                     | N/A<br>AD4500                 | W2372 (Active)<br>W4472BB Dispersive                           |   |
|                           | HV 20                        |                               | W2672 (Red Active)<br>W2772 (Blk Active)<br>W2372 (Dispersive) |   |
|                           | IF 200                       | AD8400                        | W8072  |   |
|                           | IF 2000                      | AD8500                        | W6572  |   |
|                           | III, IIIG                    | AD8500                        | W6472  |   |
|                           | RM6                          | AD8500 with<br>or AD9500 with | W6572<br>W7572   |   |
|                           | Theramini 2, 3C, 3P          | AD8500 with<br>or AD9500 with | W6572<br>W7572   | W4772 (Ban.-Pin) ,<br>W4872 (Ban.-Ban.) |
|                           | Theratach 4.7, 7.7           | AD8500 with<br>or AD9500 with | W6572<br>W7572   | W4772 (Ban.-Pin) ,<br>W4872 (Ban.-Ban.) |
| VI                        |                              | W2972RB                       |  |   |
| Scrip                     | Superstim MF Model 4000      | AD8500                        | W6572  | W7272IVRR                               |
| Sterne                    | Acutrain V                   | AD8500                        | W6572<br>W5672 (Cdn. only)                                     | W7272IVRR                               |
|                           | Myotrain V                   | AD8500                        | W6572<br>W5672 (Cdn. only)                                     | W7272IVRR                               |
|                           | NeuroTrain III               |                               |  | W7272                                   |
| US Medical                | Model 450                    |                               | W2372  |   |
| Uniphy                    |                              |                               |  | W4872                                   |
| Vectrodyne                |                              | AD8400                        | W8072  |   |
| Verimed                   | Biofeedback                  |                               |  | W2172                                   |
| Williams                  | Intertron 4000 etc.          | AD6500 with<br>or AD7500 with | W6572<br>W7572   | W8172                                   |
|                           | ZC2, ZS2, ZS4                | AD6500 with<br>or AD7500 with | W6572<br>W7572   | W8172                                   |